



American-Made Quality Meets State-of-the-Art Design

2022

American Made

All chairs, lights, delivery systems and supporting equipment are engineered and manufactured under our own roof in Hillsboro, Oregon.

Customizable

Our products are designed to complement any office design, space constrictions, functional requirements and price point.

World-Class Service

We give you access to a team of technicians and the manufacturing team itself for fast, efficient service and peace of mind.

EARN'S ITS LIVING DOUBLES



AS A DENTAL CHAIR. AS A PATIENT PLEASER.

engle
dental systems





For the Past 50 Years,

Engle has built its reputation on quality, ingenuity, and integrity. With industry-leading engineering and customer service that dentists rely on, Engle products are American-made, designed and manufactured in Hillsboro, Oregon.

Engle Dental is the only company that makes the fully programmable, 10.25" forward/backward hydraulic Traverse Chair™. The traverse feature allows dentists to position themselves behind the patient comfortably without sacrificing space. Since the seat moves forward as the back of the chair reclines, dentists can stay put and let the patient come to them. The Traverse Chair™ keeps the patient's head/oral cavity in the same vertical plane in the operator (front to back) regardless of the chair position.

Engineered for ergonomic efficiency, the Engle Traverse Chair™ allows for unparalleled dental team access and comfort. With an anatomically correct back pivot point and synchronized back and toe board recline, the Traverse Chair™ has undergone extensive comfort testing to provide maximum comfort in procedures lasting 45 minutes or more.

With options like heat and massage, lumbar support, a double articulating headrest, and Ultraleather® upholstery, dentists can focus on patient care knowing that Engle has created a level of comfort that helps patients relax.

Table of Contents

PACKAGES.....	10
• 12 O’Clock Delivery System	10
• 360 Packages	12
- Traverse	
- Standard	
• 320 Packages.....	14
- Traverse	
- Standard	
• 310 Packages.....	16
- Traverse	
- Standard	
LIGHTS.....	18
- LED Operatory Lights	
- Monitor Assembly	
CHAIRS.....	20
UPHOLSTERY.....	22
DELIVERY SYSTEMS.....	24
- Wall Mounted Delivery	
- Hygiene Delivery	
- Mobile Delivery	
- Assistant Packages	
STOOLS.....	34
TECHNICAL DIMENSIONS.....	36

Contact Info

TECHNICAL SUPPORT

TOLL FREE: 1-800-547-1906

support@engledental.com

Available: 6:00am-3:30pm PST Monday-Friday

CUSTOMER SERVICE

TOLL FREE: 1-800-547-1906

info@engledental.com

SHIPPING ADDRESS

7205 NE Evergreen Pkwy, Suite 100
Hillsboro, OR 97124

SHIPPING INFORMATION

The authorized dealer is responsible for providing a correct shipping address and details for delivery. Products are shipped F.O.B. Hillsboro, Oregon, USA. All items should be inspected upon receipt. Any shipping damage must be filed by the receiver with the carrier.

RETURN POLICY

No products will be accepted for return without prior manufacturer authorization via RMA number. RMA can be generated from Engle's technical support team. RMA number must be marked on the outside of boxes being returned, otherwise they will be refused and returned. All returned goods are subject to a 20% restocking fee. Returned goods must be received within 30 days of RMA number. All shipping charges for return of products are the sole responsibility of the dealer.

CANCELLATION POLICY

Changes or cancellations must be received by Engle Dental within 72 hours after order acknowledgment is received. Any changes or cancellation beyond that time period will be subject to additional charges or delays to order. Contact Engle Customer Service for additional questions to specific orders. For products with custom colors, including but not limited to, dental chair upholstery, stools, dental furniture, and table tops, no return will be accepted and no credits will be available. Equipment and Parts are sold exclusively through authorized Engle Dealers. Only authorized Engle dealers are permitted to sell, service Engle Equipment and Parts. Contact Engle Customer Service for a list of Authorized Equipment Dealers.

WARRANTY

All Products except stools, sold to and installed by authorized Engle Dental Systems dealers are guaranteed against defects in materials and workmanship for a period of five (5) years after date of purchase. Stools are guaranteed against defects in materials and workmanship for a period of two (2) years after date of purchase. Upholstery is guaranteed for one (1) year after date of purchase. Handpiece illumination systems are guaranteed for six (6) months after date of purchase. This warranty does not cover light bulbs of any kind. Faulty parts as determined by dealer service representative, will be repaired or replaced; no claim for labor or consequential damages will be allowed. Invoices for replacement parts may be subject to adjustment after the manufacturer has inspected the alleged defective parts. The warranty does not cover normal wear, stains, cuts, or scratches of upholstery surface finishes, or parts sold to O.E.M. customers. This warranty does not cover damage resulting from the use of cleaning, disinfecting, or sterilization chemicals and processes. The warranty does not cover damage from improper installation or maintenance, accident, misuse, or accidents in shipping and handling. In the event of shipping damage, the purchaser must initiate all claims against the freight carrier in writing within twenty-four (24) hours after the damaged items are received. The claim is the responsibility of the customer.



Ordering Information

ORDERING INFORMATION

When ordering, please include the following information: Engle is sold exclusively through authorized Engle Dealers. To find an authorized Engle equipment dealer, call 1-800-547-1906 or email info@engledental.com

- Model and/or Part Number
- Description/Color
- Quantity
- Date Required
- Purchase Order Number
- Shipping Instructions

All prices quoted herein are in US dollars and are subject to change without notice.
All prices listed in this catalog are suggested retail prices.

CUSTOMS

When ordering, you are considered the importer of record and must comply with all laws and regulations of the country in which you are receiving the goods. Contact your local customs office for further information, as customs policies vary from country to country.

RETURNED GOODS

- Before returning the product, call Customer Service for a Return Authorization Number.
Please have the following information when you call:
 - Customer number
 - Part numbers and quantities for all return items
 - The corresponding invoice or order number
- Carefully package all return items to prevent shipping damage.
- Write the Return Authorization Number on the outside of the package and reference the number on any accompanying items.
- Customized products and upholstery returns will not be accepted.
- 20% Restocking Fee


TERMS

Net 30 days on approved accounts; Visa, Mastercard or American Express can be used until credit is approved.

WARNING

Only authorized dealers should attempt to service.
Use of other than authorized dealers will void the warranty

WHAT IS TRAVERSE?



Engle is the only company that makes the "fully programmable 10.25" forward/backward hydraulic Traverse Chair™. The traverse feature allows dentists to position themselves behind the patient comfortably without sacrificing space. Since the seat moves forward as the back of the chair reclines, you can stay put and let the patient come to you.

- FLEXIBLE FOR ANY SIZE OPERATORY
- ENHANCED PATIENT COMFORT FOR LONG PROCEDURES
- ERGONOMICALLY SUITED FOR ANY HEIGHT OF DOCTOR/ASSISTANT
- IDEALLY DESIGNED TO COMBINE WITH REAR DELIVERY FOR 12 O'CLOCK DENTISTRY




DOCTOR'S & ASSISTANT'S WITH WORK SURFACE



(cabinet is not included in pricing)


PO71438
\$6,909.09



PO71450
\$6,909.09

Part Number	Equipment Included	List Price
PO71438	Cabinet Top Mount 12 O'Clock Rear Delivery	\$6,909.09
PO71450	Cabinet Side Mount 12 O'Clock Rear Delivery	\$6,909.09

CHOOSE SYSTEM

 **WARNING:** This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

 **WARNING:** This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.





A REAL UN**CROWD**ED PLEASER

ENGLE REAR DELIVERY + ENGLE TRAVERSE CHAIR™

Enjoy **UNPARALLELED ERGONOMICS** and **10.25 EXTRA INCHES** of space behind the patient with the Engle Traverse Chair and the New Engle Rear Delivery.

For more information, contact 800.547.1906

AVAILABLE AS TRAVERSE OR STANDARD PACKAGE



ASSISTANT'S

PO22275
\$1,809.09

DELIVERY UNIT

PO70931
\$6,000.00



#P096265

- Shown w/ Optional:
- Assistant's Touch Pad
 - 2nd H.V.E
 - Bien Air Handpiece
 - Chair Foot Control
 - Metal Tray
 - Twin Stitch Upholstery
 - 360 Cuspidor pictured above

Part Number	Equipment Included	List Price
CHOOSE CHAIR		
PO96540	Standard Chair	\$8,909.09
PO96540T	Traverse Chair	\$11,272.73
PO22275	Assistants	\$1,809.09
PO70931	Delivery Unit	\$6,000.00

 **WARNING:** This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

 **WARNING:** This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

CHAIR FEATURES – 5 YEAR WARRANTY

- Std. or Traverse Model
- Traverse features 10.25" Movement
- Synchronized back and seat recline
- Dual articulating headrest
- 6 Presets on Touchpad with Entry & Exit
- 60 Degree swivel
- Wideback or Narrowback w/Slings
- Swing down armrest
- Std. Seamless ColorGuard Upholstery
- 120V Chair (220V Upon Request)

LEFT/RIGHT CONVERTIBLE DELIVERY SYSTEM FEATURES – 5 YEAR WARRANTY

- E300 Control Head w/ 3 Auto HP Holders
- Junction Box w/ Air & Water Regulator
- Flex arm with Air Brake
- Master On/Off Switch
- Asepsis HP Tubing
- 2L Water Bottle System w/QD
- 3-Way Air-Water Syringe
- City/Water Bottle Switch
- Handpiece Flush/Purge System
- Quick Connect Water Aux
- Adjustable Swivel Tray
- Individual Water Coolant Controls
- Fully integrate Electric Handpieces, Intraoral Cameras, Scalers, Digital Sensors
- Wet/Dry Foot Control
- 2 FT Umbilical


ASSISTANT FEATURES – 5 YEAR WARRANTY


- Left/Right Convertible
- 3-Way Air-Water Syringe
- H.V.E., Saliva Ejector, Solids Collector
- Upgrade w/ Porcelain Cuspidor #P069265

LED LIGHT – 5 YEAR WARRANTY

- 3-Axis Rotation
- 3 Levels of Intensity
- LUX (15,000 – 25,000 – 35,000)
- Built-in Curing Filter
- 5200K Color Temperature

Part Number	Part Description	List Price	Part Number	Part Description	List Price
TRAVERSE PKG OPTIONS			STANDARD PKG OPTIONS		
P096731	Touchpad - Left Side Chair	\$590.91	P096662	Touchpad - Left Side Chair	\$536.36
P096730	Touchpad - Right Side Chair	\$590.91	P096660	Touchpad - Right Side Chair	\$536.36
96-735-01	Touchpad (Control Head Mounted-E300)	\$590.91	96-530-01	Touchpad (Control Head Mounted-E300)	\$536.36
96-735-07	Touchpad (Cabinet Mount)	\$590.91	96-530-07	Touchpad (Cabinet Mount)	\$536.36
96-735-08	Touchpad (360 Asst, P022275)	\$590.91	96-530-08	Touchpad (360 Asst, P022275)	\$536.36
P096780T	Chair Foot Control	\$772.73	P096780	Chair Foot Control (Standard)	\$700.00
P099031	Tech Connect USB Upgrade (2 USB Ports & Sensor Holder)	\$681.82	P099031	Tech Connect USB Upgrade (2 USB Ports & Sensor Holder)	\$681.82
P099033	Audio Jack Integrated into Chair	\$272.73	P096456	I.V. Arm Pad	\$536.36
P096456	I.V. Arm Pad	\$536.36	P099018	Add 4th Handpiece Position	\$727.27
P099018	Add 4th Handpiece Position	\$727.27	H014356	6 PIN Handpiece Illumination	\$545.45
H014356	6 PIN Handpiece Illumination	\$545.45	H014308	5 PIN Handpiece Illumination	\$545.45
H014308	5 PIN Handpiece Illumination	\$545.45	P070370M	Stainless Steel Metal Tray	\$154.55
P070370M	Stainless Steel Metal Tray	\$154.55	P022755	2nd H.V.E	\$254.55
P022755	2nd H.V.E	\$254.55			

 **WARNING:** This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

 **WARNING:** This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

SIDE DELIVERY OPERATORY PACKAGE



DELIVERY UNIT
W/ ASSISTANTS

PO70925
\$7,818.18

Part Number	Equipment Included	List Price
CHOOSE CHAIR		
PO96540	Standard Chair	\$8,909.09
PO95450T	Traverse Chair	\$11,272.73
PO70925	Delivery Unit	\$7,818.18

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.



CHAIR FEATURES – 5 YEAR WARRANTY

- Std. or Traverse Model
- Traverse features 10.25" Movement
- Synchronized back and seat recline
- Dual articulating headrest
- 6 Presets on Touchpad with Entry & Exit
- 60 Degree swivel
- Wideback or Narrowback w/Slings
- Swing down armrest
- Std. Seamless ColorGuard Upholstery
- 120V Chair (220V Upon Request)


DELIVERY SYSTEM FEATURES – 5 YEAR WARRANTY


- E300 Control Head w/ 3 Auto HP Holders
- Utility Center w/quick connect Water Aux
- Flex arm with Air Brake
- Master On/Off Switch
- Asepsis HP Tubing
- 2L Water Bottle System w/QD
- 3-Way Air-Water Syringe
- City/Water Bottle Switch
- Handpiece Flush/Purge System
- Wet/Dry Foot Control
- Adjustable Swivel Tray
- Individual Water Coolant Controls
- 5 FT Umbilical
- Junction Box w/ Air & Water Regulator
- H.V.E, Saliva Ejector, Solids Collector
- Telescopic Assistants Arm w/4 Pos. Holder
- Fully integrate Electric Handpieces, Intraoral Cameras, Scalers, Digital Sensors
- Assistant's 3-way Air-Water Syringe

LED LIGHT – 5 YEAR WARRANTY

- 3-Axis Rotation
- 3 Levels of Intensity
- LUX (15,000 – 25,000 – 35,000)
- Built-in Curing Filter
- 5200K Color Temperature

Part Number	Part Description	List Price	Part Number	Part Description	List Price
TRAVERSE PKG OPTIONS			STANDARD PKG OPTIONS		
P096731	Additional Touchpad - Left Side Chair (Traverse)	\$590.91	P096662	Additional Touchpad - Left Side Chair (Non-Traverse)	\$536.36
P096730	Additional Touchpad - Right Side Chair (Traverse)	\$590.91	P096660	Additional Touchpad - Right Side Chair (Non-Traverse)	\$536.36
96-735-01	Additional Touchpad (Control Head Mounted-E300)	\$590.91	96-530-01	Additional Touchpad (Control Head Mounted-E300)	\$536.36
96-735-07	Additional Touchpad (Cabinet Mount)	\$590.91	96-530-07	Additional Touchpad (Cabinet Mount)	\$536.36
96-735-08	Additional Touchpad (360 Asst, P022275)	\$590.91	96-530-08	Additional Touchpad (360 Asst, P022275)	\$536.36
P096780T	Chair Foot Control (Traverse)	\$772.73	P096780T	Chair Foot Control (Standard)	\$700.00
P099031	Tech Connect USB Upgrade (2 USB Ports & Sensor Holder)	\$681.82	P099031	Tech Connect USB Upgrade (2 USB Ports & Sensor Holder)	\$681.82
P099033	Audio Jack Integrated into Chair	\$272.73	P096456	I.V. Arm Pad	\$536.36
P096456	I.V. Arm Pad	\$536.36	P099018	Add 4th Handpiece Position	\$727.27
P099018	Add 4th Handpiece Position	\$727.27	H014356	6 PIN Handpiece Illumination	\$545.45
H014356	6 PIN Handpiece Illumination	\$545.45	H014308	5 PIN Handpiece Illumination	\$545.45
H014308	5 PIN Handpiece Illumination	\$545.45	P070370M	Stainless Steel Metal Tray	\$154.55
P070370M	Stainless Steel Metal Tray	\$154.55	P022755	2nd H.V.E	\$254.55
P022755	2nd H.V.E	\$254.55			

 **WARNING:** This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

 **WARNING:** This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

OVER THE PATIENT WITH CUSPIDOR OPERATORY PACKAGE



Shown w/ Optional:
• 2nd H.V.E

**DELIVERY
UNIT**

PO70920
\$9,081.82



Shown w/ Optional:
• Metal Tray P070370M



Shown w/ Optional:
• H.V.E

Part Number	Equipment Included	List Price
CHOOSE CHAIR		
PO96540	Standard Chair	\$8,909.09
PO95450T	Traverse Chair	\$11,272.73
PO70920	Delivery Unit	\$9,081.82

⚠ WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

⚠ WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.



CHAIR FEATURES – 5 YEAR WARRANTY

- Std. or Traverse Model
- Traverse features 10.25" Movement
- Synchronized back and seat recline
- Dual articulating headrest
- 6 Presets on Touchpad with Entry & Exit
- 60 Degree swivel
- Wideback or Narrowback w/Slings
- Swing down armrest
- Std. Seamless ColorGuard Upholstery
- 120V Chair (220V Upon Request)

DELIVERY SYSTEM FEATURES – 5 YEAR WARRANTY

- E300 Control Head w/ 3 Auto HP Holders
- Utility Center w/quick connect Water Aux
- Porcelain Cuspidor w/ Cup Fill & Bowl Rinse
- Master On/Off Switch
- Flex arm with Air Brake
- 2L Water Bottle System w/QD
- Asepsis HP Tubing
- City/Water Bottle Switch
- 3-Way Air-Water Syringe
- Wet/Dry Foot Control
- Handpiece Flush/Purge System
- Individual Water Coolant Controls
- Adjustable Swivel Tray
- Junction Box w/ Air & Water Regulator
- 5 FT Umbilical
- Telescopic Assistants Arm w/4 Pos. Holder
- H.V.E, Saliva Ejector, Solids Collector
- Assistant's 3-way Air-Water Syringe
- Fully integrate Electric Handpieces, Intraoral Cameras, Scalers, Digital Sensors

LED LIGHT – 5 YEAR WARRANTY

- 3-Axis Rotation
- 3 Levels of Intensity LUX (15,000 – 25,000 – 35,000)
- Built-in Curing Filter
- 5200K Color Temperature

Part Number	Part Description	List Price	Part Number	Part Description	List Price
TRAVERSE PKG OPTIONS			STANDARD PKG OPTIONS		
P096731	Additional Touchpad - Left Side Chair (Traverse)	\$590.91	P096662	Additional Touchpad - Left Side Chair (Non-Traverse)	\$536.36
P096730	Additional Touchpad - Right Side Chair (Traverse)	\$590.91	P096660	Additional Touchpad - Right Side Chair (Non-Traverse)	\$536.36
96-735-01	Additional Touchpad (Control Head Mounted-E300)	\$590.91	96-530-01	Additional Touchpad (Control Head Mounted-E300)	\$536.36
96-735-07	Additional Touchpad (Cabinet Mount)	\$590.91	96-530-07	Additional Touchpad (Cabinet Mount)	\$536.36
96-735-08	Additional Touchpad (360 Asst, P022275)	\$590.91	96-530-08	Additional Touchpad (360 Asst, P022275)	\$536.36
P096780T	Chair Foot Control (Traverse)	\$772.73	P096780T	Chair Foot Control (Standard)	\$700.00
P099031	Tech Connect USB Upgrade (2 USB Ports & Sensor Holder)	\$681.82	P099031	Tech Connect USB Upgrade (2 USB Ports & Sensor Holder)	\$681.82
P099033	Audio Jack Integrated into Chair	\$272.73	P096456	I.V. Arm Pad	\$536.36
P096456	I.V. Arm Pad	\$536.36	P099018	Add 4th Handpiece Position	\$727.27
P099018	Add 4th Handpiece Position	\$727.27	H014356	6 PIN Handpiece Illumination	\$545.45
H014356	6 PIN Handpiece Illumination	\$545.45	H014308	5 PIN Handpiece Illumination	\$545.45
H014308	5 PIN Handpiece Illumination	\$545.45	P070370M	Stainless Steel Metal Tray	\$154.55
P070370M	Stainless Steel Metal Tray	\$154.55	P022755	2nd H.V.E	\$254.55
P022755	2nd H.V.E	\$254.55			

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.



STANDARD FEATURES

- 3 Axis Rotation Light Head
- Focal Point: 3"x 6" @ 27"
- Low: 15,000 LUX
Medium: 24,000 LUX
High: 35,000 LUX
- Sliding Curing Filter
- No Infrared Light
- 5200K Color Temperature
- 10+ Years of LED Life
- Low Voltage

Part Number	Part Description	List Price
LED LIGHTS		
P024730	360 MAIA LED LIGHT (RADIUS), Must Be Ordered With Engle Chair	\$3,600.00
P024720	320/310 MAIA LED Light, Utility Center Mounted	\$3,600.00
P024710	MAIA LED Post Mounted Light (Need Chair Bracket & Post)	\$3,600.00
P024145	32" Engle Light Post (P024710)	\$245.45
P024140	40" Engle Light Post (P024710)	\$245.45
P024681	Antipolymerization Filter (Curing Filter)	\$163.64

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.




STANDARD FEATURES

- HDMI Cable Factory Installed
- Mounts to Engle Chair or Utility Center
- VESA Mount Compatible, size 100mm and 75mm
- Articulating Monitor Mount with Quick Click Positioning
- Monitor Recommendations:
 - 27" Max Monitor Size
 - DC Power Only

** P024620:
Monitor Not Included



P099032



P024610

Part Number

Part Description

List Price

MONITOR MOUNTS

P024610
P099032

360 Monitor Mount (Radius)
320/310 Utility Center Mounted Monitor

\$2,181.82
\$1,272.73

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

ENGLE 300 & ENGLE 300 TRAVERSE




525LB LIFT CAPACITY –
ALL HYDRAULIC MOVEMENT



P096540 - Standard Chair
P096540T - Traverse Chair

10.25" TRAVERSE FORWARD/BACK

Part Number	Part Description	List Price
CHOOSE CHAIR		
P096540T	Engle Traverse Hydraulic Chair	\$11,272.73
P096540	Engle Standard Hydraulic Chair	\$8,909.09

 **WARNING:** This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

 **WARNING:** This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.



STANDARD FEATURES


STANDARD CHAIR


- Thin Back w/ Lumbar Support
- Wide Back OR Narrowback w/ Slings
- Seamless ColorGuard Upholstery
- Anatomically Correct Back Pivot Point
- 6 Programmable Positions/4 Manual
- Touchpad (Specify Location - Lt or Rt Side)
- Synchronized Back & Toe Board Recline
- Smooth, Quiet All Hydraulic System
- 60° Seat Rotation
- Quick Release Double Articulating Headrest
- Left/Right Swing Down Upholstered Arms
- Programmable Exit

TRAVERSE CHAIR

- Thin Back w/ Lumbar Support
- Wide Back OR Narrowback w/ Slings
- Seamless ColorGuard Upholstery
- Anatomically Correct Back Pivot Point
- 6 Programmable Positions/6 Manual
- Touchpad (Specify Location - Lt or Rt Side)
- Synchronized Back & Toe Board Recline
- Smooth, Quiet All Hydraulic System
- 60° Seat Rotation
- Quick Release Double Articulating Headrest
- Left/Right Swing Down Upholstered Arms
- Programmable Exit
- 10.25" Hydraulic Traverse Movement

Part Number	Part Description	List Price	Part Number	Part Description	List Price
TRAVERSE CHAIR OPTIONS			STANDARD CHAIR OPTIONS		
P096731	Additional Touchpad - Left Side Chair Mount (Traverse)	\$590.91	P096662	Additional Touchpad - Left Side Chair Mount (Non-Traverse)	\$536.36
P096730	Additional Touchpad - Right Side Chair Mount (Traverse)	\$590.91	P096660	Additional Touchpad - Right Side Chair Mount (Non-Traverse)	\$536.36
96-735-01	Additional Touchpad (Control Head Mounted-E300)	\$590.91	96-530-01	Additional Touchpad (Control Head Mounted-E300)	\$536.36
96-735-07	Additional Touchpad (Cabinet Mount)	\$590.91	96-530-07	Additional Touchpad (Cabinet Mount)	\$536.36
P096780T	Chair Foot Control (Traverse)	\$772.73	P096780	Chair Foot Control (Standard)	\$700.00
P099031	Tech Connect USB Upgrade (2 USB Ports & Sensor Holder)	\$681.82	P099031	Tech Connect USB Upgrade (2 USB Ports & Sensor Holder)	\$681.82
P099033	Audio Jack Integrated into Chair	\$272.73	P099033	Audio Jack Integrated into Chair	\$272.73
P096456	I.V. Arm Pad	\$536.36	P096456	I.V. Arm Pad	\$536.36
P099018	Add 4th Handpiece Position	\$727.27	P099018	Add 4th Handpiece Position	\$727.27
H014356	6 PIN Handpiece Illumination	\$545.45	H014356	6 PIN Handpiece Illumination	\$545.45
H014309	5 PIN Handpiece Illumination	\$545.45	H014309	5 PIN Handpiece Illumination	\$545.45
PO70370M	Stainless Steel Metal Tray	\$154.55	PO70370M	Stainless Steel Metal Tray	\$154.55
PO22755	2nd H.V.E	\$254.55	PO22755	2nd H.V.E	\$254.55

 **WARNING:** This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

 **WARNING:** This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

STANDARD ULTRALEATHER COLOR OPTIONS

Ultra leather is the ultimate in softness and durability. Performance fabric redefined.
Engle offers 35 Standard Colors along with the full line (100+) of "Custom ULTRALEATHER" colors.

						
Buff	Cashmere	Pecan	Papyrus	Stone	Curry	Artic
						
Dove Grey	Granite	Schooner	Raven Wing	Blue Bird	Cyan	Lagoon
						
Blueberry	Diplomat Blue	Violet	Baltic	Plum	Nile	Admiral
						
Aubergine	Sorbet	Chianti	Sassy	Apricot	Campfire	Brick
						
Eucalyptus	Shetland	Parrot	Oz	Taupe	Fudge	

STANDARD COLORGUARD OPTIONS

ColorGuard features a durable, 36 oz. vinyl construction that is resistant to microbial and bacterial attack. Featuring the industry-leading PreFixx protective finish.
Engle offers 10 Standard Colors along with the full line (70+) of "Custom ColorGuard" colors.

						
Jade	Beige	Black	Coffee Cream	Colonial Blue	Gris	Lite Gray
			<div>Request Free Upholstery Samples Online at www.engledental.com or email info@engledental.com Please note printed catalog and computer images of color shades are approximate only, actual material swatch will vary in shade/tone from colors shown.</div>			
Midnight Blue	Natural	Black Plum				



Twin Switch



Summit



Seamless

Part Number**Part Description****List Price****CHOOSE UPHOLSTERY: OPTION OF WIDEBACK OR NARROWBACK W/ SLINGS***

*NARROWBACK AVAILABLE WITHOUT SLINGS UPON REQUEST

P096683	Wideback - Seamless ColorGuard Upholstery	\$1,090.91
P096685	Wideback - Upgrade to Seamless Ultraleather Upholstery	\$1,627.27
P096687	Wideback - Upgrade to Twin Stitch Ultraleather Upholstery	\$1,809.09
P096824	Wideback - Upgrade to Summit Ultraleather Upholstery	\$1,772.73
P096685HM	Wideback - Upgrade to Seamless Heat & Massage Ultraleather Upholstery	\$2,536.36
P096687HM	Wideback - Upgrade to Twin Stitch Heat & Massage Ultraleather Upholstery	\$2,900.00
P096825HM	Wideback - Upgrade to Summit Heat & Massage Ultraleather Upholstery	\$2,863.64
P096684	Narrowback w/Slings - Seamless ColorGuard Upholstery	\$1,090.91
P096686	Narrowback w/Slings - Upgrade to Seamless Ultraleather Upholstery	\$1,627.27
P096688	Narrowback w/Slings - Upgrade to Twin Stitch Ultraleather Upholstery	\$1,809.09
P096826	Narrowback w/Slings - Upgrade to Summit Ultraleather Upholstery	\$1,772.73
P096686HM	Narrowback w/Slings - Upgrade to Seamless Heat & Massage Ultraleather Upholstery	\$2,536.36
P096688HM	Narrowback w/Slings - Upgrade to Twin Stitch Heat & Massage Ultraleather Upholstery	\$2,900.00
P096827HM	Narrowback w/Slings - Upgrade to Summit Heat & Massage Ultraleather Upholstery	\$2,863.64
EMBR	Custom Embroidery - Image file Required with Order	\$1,091.91
DUAL	Dual Tone or Custom Upholstery	\$1,454.55

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.



Part Number

Part Description

List Price

Part Number

Part Description

List Price

DELIVERY SYSTEM

ADDITIONAL OPTIONS

P070936	Engle Dr's Side Delivery System, Wall Mount, E300 Control Head	\$4,909.09
P070936S	Engle Dr's Side Delivery System, Wall Mount (Short Flex Arm),	\$4,909.09
P070900	AS-1 Control Head Engle Dr's & Asst Delivery System, Wall Mount	\$3,418.18
P070902	Engle Dr's Wall Mount Delivery System, AS-1 Control Head	\$2,909.09
P029355	Engle Panel Mount Control System	\$2,272.73
P070904	Engle Assistant's Delivery System, Wall Mount	\$1,090.91

H014356	6 PIN Handpiece Illumination	\$545.45
H014308	5 PIN Handpiece Illumination	\$545.45
P070370M	Stainless Steel Tray	\$154.55
P071655	Wall Mounting Board	\$272.73
P021604	Clean Water Auxiliary System w/Regulator	\$336.36

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

STANDARD FEATURES**P070902**

- AS-1 Control Head w/ 3 Auto Handpiece Positions
- Individual Water Coolant Controls
- Asepsis Handpiece Tubing
- Wet/Dry Foot Control
- Handpiece Flush/Purge System
- 3-Way Air-Water Syringe
- Master On/Off Switch
- Small Junction box w/ Air & Water Master Valve /Filter/Regulator
- Smooth 3-Pivot Articulating Control Head Delivery
- Individual Drive Air Adjustments Per Handpiece
- Left/Right Delivery
- 4' Long Umbilical

P070904

- Asepsis Handpiece Tubing
- Solids Collector
- H.V.E.
- Saliva Ejector
- 4-Pivot Articulating Assistant's Arm w/ 4-Position Holder
- Left/Right Delivery
- 4' Long Umbilical
- Assistants 3-way Air/Water Syringe

STANDARD FEATURES**P070936**

- E300 Control Head w/ 3 Auto Handpiece Holders
- Individual Water Coolant Controls
- Master On/Off Switch
- Asepsis Handpiece Tubing
- Handpiece Flush/Purge System
- 3"-Way Air-Water Syringe

P070936S

- E300 Control Head w/ 3 Auto Handpiece Positions
- Short Flex Arm
- Flex Arm w/ Air Brake
- Adjustable Swivel Tray
- Small Junction Box w/ Air & Water Master Valve/Filter/Regulator Assembly
- 8' Long Umbilical
- Wet/Dry Foot Control

STANDARD FEATURES**P070900**

- AS-1 Control Head w/ 3 Auto Handpiece Positions
- Individual Water Coolant Controls
- Asepsis Handpiece Tubing
- Wet/Dry Foot Control
- Handpiece Flush/Purge System
- 3-Way Air-Water Syringe
- Master On/Off Switch
- Solids Collector
- H.V.E.
- Saliva Ejector
- Small Junction box w/ Air & Water Master Valve /Filter/Regulator Assembly
- 4-Pivot Articulating Assistant's
- Arm w/ 4-Position Holder
- Smooth 3-Pivot Articulating Control Head Delivery
- Individual Drive Air Adjustments Per Handpiece
- Left/Right Delivery
- 4' Long Umbilical
- Assistants' #-way Air/Water Syringe

P029355

- Panel Mount Control w/ 3 Auto Handpiece Holders
- Individual Water Coolant Controls
- Asepsis Handpiece Tubing
- Wet/Dry Foot Control
- Handpiece Flush/Purge System
- 3-Way Air-Water Syringe
- Master On/Off Switch
- Air & Water Master Valve/Filter/Regulator Assembly
- Individual Drive Air Adjustments Per Handpiece
- 6' Long Umbilical

STANDARD FEATURES

- P071964**
- Under Counter Mounted AS-1 Delivery Head
 - Left / Right Delivery
 - 3 Automatic Handpiece controls
 - Individual Water Coolant Controls
 - Individual drive air controls
 - Asepsis Handpiece Tubing
 - Wet/Dry Foot Control
 - Handpiece Flush/Purge System
 - 3-Way Air/Water Syringe
 - Master On/Off Switch
 - 4' Umbilical



Part Number	Part Description	List Price
DELIVERY SYSTEM		
P071964	Rear Delivery, Under Counter Mount - Doctor	\$3,272.73
DELIVERY SYSTEM OPTIONS		
H014356	6 Pin Handpiece Illumination	\$545.45
H014308	5 Pin Handpiece Illumination	\$545.45
PO21604	Clean Water Auxiliary W/ Regulator	\$336.36

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.





P070937



P070886



P070885

Part Number	Part Description	List Price
DELIVERY SYSTEM		
P070928	Engle 360 Delivery: Radius, Left/Right AS-1 Control Head Delivery System	\$5,454.55
P070885	Engle 1200 Delivery: Over the Patient, AS-1 Control Head Delivery System w/Cuspidor	\$8,545.45
P070886	Engle 2200 Delivery: Side Mount, AS-1 Control Head Delivery System	\$6,181.82
P070912	Engle 2200 Delivery: Chair Mount, AS-1 Control Head, without Assistant's System	\$5,454.55
P070937	Engle 360 Hygiene Delivery: Radius, Left/Right, AS-1 Control Head w/Vacuum Delivery System (Mounts to Standard Chair Only)	\$5,272.73

** P070937 NOT COMPATIBLE WITH TRAVERSE CHAIR

STANDARD FEATURES

DELIVERY

- AS-1 Control Head With 2 Automatic Handpieces Saliva Ejector And H.V.E.
 - Left/Right Configuration
 - Individual Water Coolant
 - Master On/Off Switch
 - Aspesis Handpiece Tubing
 - Handpiece Flush/Purge
- Wet/Dry Foot Control
 - 3-way Air/Water Syringe
 - Saliva Ejector
 - H.V.E.- Solids Collector
 - Flex Arm W/ Air Brakes
 - 2 Liter Water Bottle With Quick Disconnect
 - City Or Water Bottle Selector
 - Contoured Box

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

DELIVERY CARTS FOR “LEAN” AND CABINET FREE OPERATORIES

STANDARD FEATURES

PLATINUM

- E300 Control Head w/3 Auto HP Holders
- Left/Right Delivery
- StarBoard Work Surface
- 3 Pivot Articulating Control Head
- 3 Pivot Articulating Asst Arm
- Master On/Off Switch
- Asepsis HP Tubing
- Wet/Dry Foot Control
- 3 Way Air-Water Syringe
- H.V.E, Saliva Ejector, Solids Collector
- Individual Water Coolant Systems
- 2L Bottled Water System
- Quick Connect Water Aux
- 4 Storage Bins w/Lid
- Junction Box, 5ft Umbilical
- 5 Year Warranty



Part Number	Part Description	List Price	Part Number	Part Description	List Price
MOBILE DELIVERY SYSTEMS			ADDITIONAL OPTIONS		
P050300	SILVER Mobile Cart DR/ASST	\$5,818.18	P051042	Corian Solid Surface 18"x30"	\$2,180.00
P050310	GOLD Mobile Cart DR/ASST	\$6,545.45	H014356	6 Pin Handpiece Illumination	\$545.45
P050320	PLATINUM Mobile Cart DR/ASST	\$7,272.73	H014308	5 Pin Handpiece Illumination	\$545.45
P050330	PLATINUM Mobile Cart ASST Only	\$6,363.64	P022755	2nd H.V.E	\$254.55
			P070370M	Stainless Steel Tray	\$154.55

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

P050135

STANDARD FEATURES

- AS-1 Control Head
- 3-Positions Automatic Handpiece
- Large Work Surface
- Individual Water Coolant Controls
- Asepsis Handpiece Tubing
- Wet/Dry Foot Control
- Handpiece Flush/Purge System
- 3-Way Air-Water Syringe
- Master On/Off Switch
- Solids Collector
- H.V.E.
- Saliva Ejector
- 2 Liter Bottled Water System
- Small Junction Box Air & Water Master Valve/Filter/Regulator Assembly
- Quick Connect Water Auxiliary
- 3-Pivot Articulating Assistant's Arm w/ 4-Position Holder
- Smooth 3-Pivot Articulating Control Head Delivery
- Left/Right Delivery
- 4 Storage Bins w/ Lid
- 5' Long Umbilical
- Assistant's 3-Way Syringe



Part Number	Part Description	List Price	Part Number	Part Description	List Price
DELIVERY CARTS - DOCTOR'S			DELIVERY SYSTEM OPTIONS		
P050135	Delivery Cart-Doctor & Assistant w/ Table Top	\$5,272.73	H014356	6 Pin Handpiece Illumination	\$545.45
			H014308	5 Pin Handpiece Illumination	\$545.45
			P099019	Add 4th Automatic Handpiece Position	\$727.27
			P021619	Quick Release For Water Bottle	\$200.00
			P022755	2nd H.V.E.	\$254.55
			P049001	Swivel Tray	\$354.55
			P070370M	Metal Tray Upgrade	\$154.55

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

P050140



STANDARD FEATURES

- 3-Position Automatic Handpiece Positions
- Individual Water Coolant Controls
- Asepsis Handpiece Tubing
- Wet/Dry Foot Control
- Handpiece Flush/Purge System
- 3-Way Air-Water Syringe
- Master On/Off Switch
- Bottled Water System
- Small Junction Box Air & Water Master Valve/Filter/Regulator Assembly
- Quick Connect Water Auxiliary
- 5' Long Umbilical

Part Number	Part Description	List Price	Part Number	Part Description	List Price
DELIVERY CARTS - DOCTOR'S			DELIVERY SYSTEM OPTIONS		
P050140	Delivery Cart - Doctor's Only	\$3,454.55	H014356	6 Pin Handpiece Illumination	\$545.45
			H014308	5 Pin Handpiece Illumination	\$545.45
			P099019	Add 4th Automatic Handpiece Position	\$727.27

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

P050145



Shown w/ Optional:
• Metal Tray

STANDARD FEATURES

- Asepsis Handpiece Tubing
- Large Work Surface
- Master On/Off Switch
- Solids Collector
- H.V.E.
- Saliva Ejector
- Bottled Water System
- Small Junction Box Air & Water Master Valve/Filter/Regulator Assembly
- Quick Connect Water Auxiliary
- 3-Pivot Articulating Assistant's Arm w/ 4-Position Holder
- Left/Right Delivery
- 4 Storage Bins w/ Lid
- 5' Long Umbilical
- Assistant's 3-Way Syringe

Part Number	Part Description	List Price	Part Number	Part Description	List Price
DELIVERY CARTS			DELIVERY SYSTEM OPTIONS		
P050145	Delivery Cart - Assistant's w/ Table Top	\$3,454.55	P022755	2nd H.V.E. for Assistant	\$254.55
			P070370M	Metal Tray	\$154.55
			P049001	Swivel Tray	\$354.55

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

P050150



STANDARD FEATURES

- 3-Position Automatic Handpiece Holders
- Individual Water Coolant Controls
- Asepsis Handpiece Tubing
- Wet/Dry Foot Control
- Handpiece Flush/Purge System
- 3-Way Air-Water Syringe
- Master On/Off Switch
- Solids Collector
- H.V.E.
- Saliva Ejector
- Bottled Water System
- Small Junction Box Air & Water Master Valve/Filter/Regulator Assembly
- Quick Connect Water Auxiliary
- 5' Long Umbilical

Part Number	Part Description	List Price	Part Number	Part Description	List Price
DELIVERY CARTS			DELIVERY SYSTEM OPTIONS		
P050150	Delivery Cart - Hygiene	\$3,454.55	H014356	6 Pin Handpiece Illumination	\$545.45
			H014308	5 Pin Handpiece Illumination	\$545.45
			P099019	4th Automatic Handpiece Position	\$727.27

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.



P069201

P069240 -
shown with optional
P020336 - Asst Syringe
P022755 - 2nd H.V.E



22-507-02



Part Number

Part Description

List Price

ASSISTANT'S AND OPTIONS

P069201	Utility Center w/ Assistant - Chair Bracket Needed	\$3,454.55
P022755	2nd H.V.E. For Assistant	\$254.55
P024343	Engle Chair Bracket	\$454.55
P069240	Utility Center w/ E300 Cuspidor and Assistant- Chair Bracket Needed	\$4,000.00
P020336	Assistant 3-WAY SYRINGE	\$181.82
P099016	AIR QD	\$163.63
P070904	Assistant's 12 O'clock Delivery System - wall / Cabinet Mounted	\$1,090.91
P022275	Engle 360 Assistant's Vacuum Delivery System - Left/Right	\$1,809.09
P069265	Engle 360 E300 Cuspidor Delivery System - Left/Right (Replaces P022275)	\$4,000.00
22-502-02	Cabinet Vacuum Group - Top Mount	\$1,090.91
22-507-01	Cabinet Vacuum - Short Arm	\$1,090.91
22-507-02	Cabinet Vacuum - Long Arm	\$1,090.91

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

**STANDARD - P097092**

- Comfortable 4" Thick, 15" Diameter Seat Cushion
- Many Options Allow Ability to Customize
- Contour Backrest (Optional)
- Various Height Options
- Swing Arm or Ratched Body Support

**CONTOUR - P097094**

- Contoured 18" x 18" Seat Cushion
- 24" Powder-Coated Cast Aluminum Base
- Dual Adjustable Lumbar Backrest Adds Comfort
- Modern Design Fits Any Office

**ERGONOMIC - P097096**

- Accommodates Front to Back Reaching without Sacrificing Lumbar Support
- Exclusive Use of HybriGel Foam
- Independently Adjust Backrest to All Greater Range of Movement
- Provides Correct Weight Distribution

**STANDARD - P097093**

- Comfortable 4" Thick, 15" Diameter Seat Cushion
- Many Options Allow Ability to Customize
- Contour Backrest (Optional)
- Various Height Options
- Swing Arm or Ratched Body Support

**CONTOUR - P097095**

- Contoured 18" x 18" Seat Cushion
- 24" Powder-Coated Cast Aluminum Base
- Dual Adjustable Lumbar Backrest Adds Comfort
- Modern Design Fits Any Office

**ERGONOMIC - P097097**

- Accommodates Front to Back Reaching without Sacrificing Lumbar Support
- Exclusive Use of HybriGel Foam
- Independently Adjust Backrest to All Greater Range of Movement
- Provides Correct Weight Distribution

Part Number**Part Description****List Price****DOCTOR & ASSISTANT STOOLS**

P097092	Doctor Stool - Standard, Vinyl Upholstery	\$1,009.09
P097093	Assistant Stool- Standard, Vinyl Upholstery	\$1,270.91
P097094	Contour Doctor Stool - Vinyl Upholstery	\$1,227.27
P097095	Contour Assistant Stool- Vinyl Upholstery	\$1,500.00
P097096	Ergonomic Doctor Stool- Vinyl Upholstery	\$2,000.00
P097097	Ergonomic Assistant Stool - Vinyl Upholstery	\$2,909.09
P097092UL	Doctor Stool - Standard, Ultraleather Upholstery	\$1,190.91
P097093UL	Assistant Stool- Standard, Ultraleather Upholstery	\$1,452.73
P097094UL	Contour Doctor Stool - Ultraleather Upholstery	\$1,409.09
P097095UL	Contour Assistant Stool- Ultraleather Upholstery	\$1,981.82
P097096UL	Ergonomic Doctor Stool- Ultraleather Upholstery	\$2,181.82
P097097UL	Ergonomic Assistant Stool - Ultraleather Upholstery	\$3,090.91

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.




ENGINEERED TO PERFORM



SERVICE & SUPPORT

Handcrafted in Hillsboro, Oregon, Engle products are engineered and manufactured under our own roof.



CUSTOMIZABLE



With access to a team of technicians and the manufacturing team itself, Engle provides fast and efficient problem-solving.

With your choice of custom colors, fabrics, and upholstery patterns, you can design your chair to suit your style.



www.engledental.com | 800.547.1906



www.engledental.com





* INDICATES FULLY TRAVERSED DIMENSION




* INDICATES FULLY TRAVERSED DIMENSION



P070931





P070350





P024667





P024610




P022275



P069265





P024617




Ceiling Post Part Number	Ceiling Post Length	Intended Ceiling Height
24-483-17	17"	8' Ceiling
24-483-23	23"	8' 6" Ceiling
24-483-29	29"	9' Ceiling
24-483-35	35"	9' 6" Ceiling
24-483-41	41"	10' Ceiling


P070925




P070920





310 UNIT ARM IDENTICAL TO 320




P024616




P070935




P070900 - DR'S D.P. WALL MOUNT




P070902 - DR'S WALL MOUNT



P070904 - ASSISTANT'S WALL MOUNT



P071964 - CABINET MOUNTED UNIT



AS SEEN IN:

DENTAL PRODUCT SHOPPER

ENGLE DENTAL SYSTEMS

Delivering a Cutting-Edge Operatory Solution

Engle's range of dental chairs, lights, delivery systems, and supporting equipment were designed to be used together to create the complete, ergonomic dental operatory

Whether you're remodeling an operatory or building a new office, choosing the right equipment can be tricky. There are dozens of factors to consider such as how much space is available; where is the best location for instruments, lighting, and monitors; and whether you practice left-, right-, or four-handed dentistry. Plus, there's the consideration of how the setup affects patient comfort, treatment



CONNECT TO THE FUTURE

Engle always has the future of digital dentistry in mind. The **TechConnect Operatory Package** is a first-of-its-kind system that powers a practice's existing digital technology through prewired HDMI cables, 2 USB ports for internal and external use, and a digital sensor holder. This offers clinicians simple chairside access to tools and technology like digital sensors, handpieces, intraoral cameras, and monitors.

time, workflow, and the overall team experience.


No matter what space and conditions your practice has to work with, Engle Dental Systems can provide customized, quality operatory equipment that meets your unique needs. For the past 50 years, Engle has developed dental chairs, delivery units, and other equipment with a focus on state-of-the-art quality and ergonomics. All equipment is designed and manufactured at its Hillsboro, OR, headquarters, and quality control is handled in-house.

Optimize the Operatory


Remodeling or designing a new room begins with the dental chair. Engle's fully programmable, 10" forward/backward Traverse Chair allows clinicians to position themselves behind the patient comfortably without sacrificing space. Built with advanced hydraulic controls, the seat moves forward as the back of the chair reclines so dentists can let the chair position itself. "The chair's ability to traverse was huge for our practice," said hygienist and *DPS* evaluator Debbie Alo, who found the Traverse Chair to be an ideal solution for her small hygiene room.



In addition to the dentist's choice of color, fabric, and upholstery pattern, the Standard and Traverse chair models feature lumbar support, a double articulating headrest, and other features to help patients relax. To complement the chair, Engle carries four operatory LED light types and a wide variety of delivery units and systems, with options for side delivery, wall-mounted, and carts. Engle customers also receive access to a world-class team of technicians and manufacturers for fast and efficient problem-solving.



PLUG INTO THE FUTURE



Digital sensor holder
on the control head

Two USB ports
pre-fitted on the
control head

Monitor mounts
pre-wired with
HDMI cables

Headphone jack
for patient
entertainment



www.engledental.com

Handcrafted in Oregon



USA

**Engle Dental Systems LLC
7205 NE Evergreen Parkway, Suite 100
Hillsboro, Oregon 97124
Tel: (800) 547-1906 • Fax: (503) 357-7286
info@engledental.com**

**Prices subject to change without notice.
P242000 6.0**