Polaroid



Polaroid Ray Digital X-Ray Sensor

Available in sizes 1.0, 1.5 and 2.0



5 - Year warranty, Free integration, Free Lifetime Support, No fees

Lowest Operational Cost over 5 years-No Fees = ROI

Minimal Purchase Investment
Free Installation and training
Free Lifetime remote support-Low cost to own
5-Yr Manufacturers Warranty-Built-in Durability
Optional 5-Yr Unconitional Accidental Breakage Policy

Integration with most Modern Imaging Softwares

Reinforced USB connector for easy direct plug in If needed, Polaroid DDS Imaging bridges to all Practice Management programs. Polaroid DDS Imaging works with ALL other Sensors, Scanners, and Pan/Cephs to improve workflow and reduce operational cost and staff learning.

Solid Construction for Durability Backed up by a 5-Year Warranty

Built to last with Internal & External impact cushions
CsI (Caesium Iodide) scintillator for the best images
Reinforced Strain relief to prevent cable damage
Reinforced Kevlar cable and USB connector
CMOS sensor for extended durability
IP-68 Waterproof rated

Superb Image Quality

High Tissue Contrast resolution

High Spacial resolution from quality CMOS sensor

High Resolution of >25lpm theoretical for sharp, defined lines

Built-In Electronic Noise and Artifact Reduction

Increased Patient Experience

Anatomical beveled corners for easy placement Rounded edges for improved patient experience Only 4.8 mm thick, Ray is the thinest sensor



Built-in Durability



Patient Safety

low energy needed for sensor activation
helps to work within the A.L.A.R.A. principles
Large Range of Clinically Acceptable Exposure levels
= Less Retakes

Choice of Sizes

All Three sizes- #1/ #1.5/ #2 are available
Our 1.5 size (Universal) helps take FMX's in half the time,
increases production, reduces cost, and improves patient comfort

Accessories

Custom molded Perfect Fit Bite Blocks
Polaroid Perfect Fit Alignment Ring and Rod kits
Perfect Fit Disposable sleeves (500Pk)

Size 2.0 dimensions	Approx 31.5 x 42.8 x 4.8 mm	
Size 2.0 active area	Approx 26 x 36 mm	
Size 1.5 dimensions	Approx 29.2 x 39.5 x 4.8 mm	
Size 1.5 active area	Approx 24 x 33 mm	
Size 1.0 dimensions	Approx 25.4 x 36.8 x 4.8 mm	
Size 1.0 active area	Approx 20 x 30 mm	
Sensor/lens type	CMOS/liquid lens Csl	
Pixel pitch	Full=14.8 Binning Mode =29.6	
A/D conversion	14 bit analog to 12 bit digital	
Interface	Direct USB	



Waterproof	Rated IP68	
USB Cable	Approx 3M (9 8ft)	
Polaroid Ray Size #1	P01-410	
Polaroid Ray Size #1.5	P01-415	181
Polaroid Ray Size #2	P01-420	
Perfect Fit #1 Alignmer	nt Ring & Rod kit w/4 Bite blocks	P01-501
Perfect Fit #1.5 Alignmen	nt Ring & Rod kit w/4 Bite blocks	P01-515
Perfect Fit #2 Alignmen	nt Ring & Rodd kit w/4 Bite blocks	P01-502
5 - Yr Accidental Damage Co	ver (Purchase with-in 30 days of sale)	P01-603

Contact your dealer for additional Polaroid imaging supplies

Dealer Contact info

Contact your local dealer for a free demo





